



# JD<sup>TM</sup> COVID-19 qRT-PCR Kit



# JD COVID-19 qRT-PCR Kit

JD COVID-19 qRT-PCR Kit is used for *in vitro* diagnostics of novel coronavirus (SARS-CoV-2) in human oropharyngeal and nasopharyngeal swab and sputum samples. JD COVID-19 qRT-PCR Kit uses Real-time RT-PCR. SARS-CoV-2 is first transcribed into cDNA by reverse transcriptase, and then cDNA is used as a template for PCR amplification. Primer and fluorescent probe were set to the N gene site in reference to US FDA's "FAQs on Diagnostic Testing for SARS-CoV-2." JD COVID-19 qRT-PCR Kit can determine positive results through reaction of increased amount of fluorescent Taqman probe contained in the reagent, in real time.

<b>Product Name</b>	JD COVID-19 qRT-PCR Kit
<b>Packaging Unit</b>	96 tests/kit, box
<b>Target Gene</b>	N gene (SARS-CoV-2)
<b>Running Time</b>	1 hour 40 minutes
<b>Storage</b>	-20°C or below
<b>Shelf-life</b>	1 year from manufacturing date
<b>Clinical Trial Results *</b>	Positive Predictive Value (PPV) = 100% (95% CI 88.78-100) Negative Predictive Value (NPV) = 100% (95% CI 90.00-100)
<b>Manufacturing</b>	GMP for Pharmaceuticals, GMP for Medical Devices, ISO 13485:2016
<b>Highly Recommended Machine</b>	Bio-Rad CFX96, Thermo Fisher Scientific AB7500, Roche LC96




\* We performed clinical trial based on 31 positive and 35 negative clinical samples in accordance with FDA's Emergency Use Authorization clinical study guidelines





# JD COVID-19 qRT-PCR Kit

## Target Gene



### Target : N gene

US CDC recommends N gene of SARS-CoV-2 for Real-Time PCR



Amplifying the target gene regardless of mutation





# JD COVID-19 qRT-PCR Kit

## Sensitivity

Diluted SARS-CoV-2 was detected with high sensitivity, with linearity.

LOD: 5 copies/μl



• Rxn cycle  
 50°C 30min  
 95°C 15min  
 95°C 5s  
 55°C 30s  
 4°C 10s  
 ↻ 40cycle

Lane	1	2	3	4	5	6	7	8	9
Name	Standard RNA	Standard RNA	Standard RNA	Standard RNA	Standard RNA	Standard RNA	Standard RNA	Standard RNA	Negative Control
Conc. (Copies/μl)	10 <sup>7</sup>	10 <sup>6</sup>	10 <sup>5</sup>	10 <sup>4</sup>	10 <sup>3</sup>	10 <sup>2</sup>	10 <sup>1</sup>	10 <sup>0</sup>	-
Result*	P	P	P	P	P	P	P	P	N
Ct values	1	13.32	16.76	20.35	23.81	27.56	30.07	34.63	37.52
	2	12.88	16.59	20.50	23.74	27.06	30.15	32.9	37.08
	3	13.34	16.91	20.51	23.91	27.23	30.65	34.89	37.50
Mean		13.18	16.75	20.45	23.82	27.28	30.29	34.14	37.37
STD		0.26	0.16	0.09	0.09	0.25	0.31	1.08	0.25
CV%		1.97	0.96	0.44	0.36	0.93	1.04	3.17	0.66


\* P: Positive, N: Negative

# JD COVID-19 qRT-PCR Kit

## Specificity

No.	Virus	Result *
1	SARS-CoV-2 (β coronavirus)	P
2	MERS-CoV (β coronavirus)	N
3	HCoV-NL63 (α coronavirus)	N
4	Human Influenza A (CA/04, H1N1)	N
5	Human Influenza H3N2	N
6	Human Influenza B (Vic)	N
7	PED(SM98) (α coronavirus)	N

\* P: Positive, N: Negative



JD COVID-19 qRT-PCR Kit  
specifically detects SARS-CoV-2.

# JD COVID-19 qRT-PCR Kit

## Instruction for Use



1. Preparation of reagent and extracted RNA
2. qRT-PCR



\* Real-time PCR Machine: CFX96 (Bio-Rad)

Step	PCR Program (20µl reaction)		
	Temp	Time	Cycle
cDNA synthesis	50°C	30 min	1 cycle
Initial inactivation	95°C	15 min	1 cycle
Denaturation	95°C	5 sec	40 cycles
Elongation	55°C	30 sec	
Hold	4°C	10 sec	1 cycle

3. Reading




4. Data Analysis and Interpretation

Item	Fluorescent		Result
	FAM(SARS-CoV-2)	HEX(Human B2M)	
SAMPLE 1	+	+/- *	Positive
SAMPLE 2	-	+	Negative
SAMPLE 3	-	-	False Negative/Retest

\* According to EUA summary (FDA), a high SARS-CoV-2 (target N gene) RNA load in the sample can lead to reduced or absent Internal Positive Control(IPC), thus making the test result positive.





\*refer to the appendix

의약품 제조 및 품질관리기준 적합판정서  
(Certificate of GMP Compliance of a Manufacturer)

- 대전지방식품의약품안전청장

의료기기 제조 및 품질관리 기준 적합인정서  
(Certificate of GMP)

- |                          |            |
|--------------------------|------------|
| 발행일자(Date of Issue)      | 2020-03-27 |
| 유효기간(Date of Expiration) | 2022-03-26 |

Enclosure  
Annoy 1

## MSCB-156

---

20

\*US FDA EUA is in progress



# JD COVID-19 qRT-PCR Kit

## F&Q

<b>MOQ</b>	1,000 Kit (96,000 tests)	
<b>MSRP</b>	Contact for pricing details	
<b>Payment Term</b>	100% T/T advanced, non-refundable	
<b>Incoterms</b>	FOB KOREA	
<b>Lead Time</b>	Within 2-week basis after payment (Firm delivery date will be notified after payment)	
<b>Certificates</b>	CFS by MFDS (KFDA) CE-IVD US FDA EUA (in process)	GMP for Pharmaceuticals GMP for Medical Devices ISO 13485:2016
<b>Storage</b>	-20 °C or below	
<b>Shelf-life</b>	1 year from manufacturing date	
<b>Weekly Manufacturing Capacity</b>	15,000 Kit	
<b>Highly Recommended PCR Machine</b>	Bio-Rad CFX96, Thermo Fisher Scientific AB7500, Roche LC96	
<b>Dimensions and Weight</b>	1kit = 70mm x 63mm x 52mm / 30g 1 container = 764mm x 590mm x 535 mm / 31,430g (±10%)	
<b>Packaging</b>	Thermal Insulation Container with Bio-Pharm Coolant (6 containers per 1,000 kits / 187,080g (±10%))	
<b>Check List for End Point</b>	Precheck for Registration or Certificate for IVD importer in destination country And confirmation of BSL(Bio Safety Level 3) lab for end-user	



# Appendix 1, KMFDS Certificate of Free Sales



Document Number : 1KHK-WJ7L-OTSO-UQRJ

Osong Health Technology Administration Complex,  
187 Osongsaengmyeong2-ro, Osong-eup, Heungdeok-gu,  
Cheongju-si, Chungcheongbuk-do, Korea, 28159  
Tel: +82-43-719-2342, Fax: +82-43-719-1000

No. of Certificate : 2020044307

Date : 2020/04/21

## Certificate of Free Sales

Exporting(certifying) country : Republic of Korea  
Importing(requesting) country :

The Ministry of Food and Drug Safety, certifies that the following firm is authorized to manufacture medical devices under the Medical Device Act and the following item(s) is(are) permitted to be freely sold in overseas markets.

Manufacturer (Registered No. : 5480)

JOONGHUN PHARMACEUTICAL Co., Ltd.

28, Oksansandan-ro, Oksan-myeon, Heungdeok-gu, Cheongju-si, Chungcheongbuk-do

Product-License No.	Classification
20-292	IVD reagents for infectious disease marker(Diagnosis of Sexually transmitted disease, Legally designated infectious pathogens other than 'high risk pathogens', infectious agents with moderate infectivity), nucleic acid test [3] Product Name : JD COVID-19 qRT-PCR Kit

\*Attached : List of Product Classification and Model

Director of High-Tech Medical Devices Division  
Department of Medical Device Evaluation  
National Institute of Food and Drug Safety  
Evaluation  
Ministry of Food and Drug Safety



\* This certificate is issued on the Internet, you can check whether to forge or modify in homepage(emed.mfds.go.kr).  
Furthermore, You can also check it by barcode exploiting document check program for scanner

## Appendix 2, CE



www.mt-procons.com

MEDICAL  
TECHNOLOGY **promedt**  
CONSULTING

# Certificate of EU product notification

Herewith we confirm that

**MT Promedt Consulting GmbH**  
**Altenhofstraße 80**  
**66386 St. Ingbert**  
**Germany**

has taken over the function of an European Authorized Representative according to the requirements of Article 10 of the IVDD 98/79/EC for

**Joonghun Pharmaceutical Co., Ltd**  
**15, Gukhoe-daero 62-gil, Yeongdeungpo-gu**  
**Seoul**  
**Republic of Korea**

MT Promedt Consulting GmbH has made the product notification at the relevant competent authority according to Article 10(3).

The in vitro diagnostic medical devices of the manufacturer, covered by the notification, are listed in Annex I of this certificate.

This certificate does not attest the conformity of the medical devices with the above mentioned directive. The conformity is stated in the respective product-related Declarations of Conformity signed under the sole responsibility of the manufacturer.

30 April 2020

Dr. Michael Rinck  
- Managing Director -

**Enclosure**  
Annex I

# Appendix 2, CE, continued



www.mt-procons.com

MEDICAL  
TECHNOLOGY **promedt**  
CONSULTING

**Joonghun Pharmaceutical Co., Ltd**

Annex I  
to "Certificate of EU Product Notification"  
(List of CE marked Products)

Page 1 / 1 of Annex I

Internal Number	Registration Number (at the German CA/ DIMDI)	Product Category (EDMS)	EDMS Code Description	Classification Annex
J00-01	DE/CA70/40838- 155295	16 90 90 01 90	Other Other Genetic Tests	III

30 April 2020

Dr. Michael Rinck  
- Managing Director -



## Appendix 3, US FDA Acknowledgment Letter



**FDA** U.S. FOOD & DRUG  
ADMINISTRATION

### Acknowledgment Letter

5/4/2020

Kyo San Ku  
AKI Health LLC  
651 N Broad St, Suite 205 #2405  
Middletown, DE 19709  
UNITED STATES

Dear Kyo San Ku:

The Center for Devices and Radiological Health (CDRH) of the Food and Drug Administration (FDA) has received your submission. This submission has been assigned the unique document control number below. All future correspondence regarding this submission should be identified prominently with the number assigned and should be submitted to the Document Control Center at the above letterhead address. Failure to do so may result in processing delays. If you believe the information identified below is incorrect, please notify the Program Operations Staff at (301) 796-5640.

Submission Number: EUA200652  
Received: 5/4/2020  
Applicant: Joonghun Pharmaceutical Co., Ltd.  
Device: JD COVID-19 qRT-PCR Kit

We will notify you when the review of this document has been completed or if any additional information is required. If you are submitting new information about a submission for which we have already made a final decision, please note that your submission will not be re-opened. For information about CDRH review regulations and policies, please refer to <http://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/default.htm>.

Sincerely yours,

Center for Devices and Radiological Health

U.S. Food & Drug Administration  
10903 New Hampshire Avenue  
Silver Spring, MD 20993  
[www.fda.gov](http://www.fda.gov)

## Appendix 4, Certificate of ISO13485

	Certificate No : IMS 6001
<b>MEDICAL DEVICE QUALITY MANAGEMENT SYSTEMS CERTIFICATE</b>	
<b>JOONGHUN PHARMACEUTICAL CO., LTD.</b>	
HEAD-OFFICE: 15, GUKHOE-DAERO 62-GIL, YEONGDEUNGPO-GU, SEOUL, 07236, KOREA FACTORY: 28, OKSANSANDAN-RO, OKSAN-MYEON, HEUNGDEOK-GU, CHEONGJU-SI, CHUNGCHEONGBUK-DO, 28101, KOREA	
This is to certify that organization above has been assessed and registered by KTR-CC for the scope of supply described as follows:	
<b>MEDICAL DEVICE QUALITY MANAGEMENT SYSTEMS</b>	
<b>ISO 13485:2016</b>	
<b>SCOPE OF SUPPLY :</b> DESIGN, DEVELOPMENT AND PRODUCTION OF GRAFT/PROSTHESIS, BIOMATERIALS AND IVD REAGENTS FOR INFECTIOUS DISEASE MARKER(JD COVID-19 qRT-PCR Kit)	
Valid Until	: FEBRUARY 27. 2022
Registered Date	: FEBRUARY 28. 2019
Issued Date	: MARCH 30. 2020
<b>KTR CERTIFICATION CENTER</b>	
#501, 646-30, Hosu-ro, Ilsandong-gu, Goyang-si, Gyeonggi-do, Korea Tel : 82-2-2093-3450/82-70-4010-4981 Web : www.ktrcc.or.kr	
 President	
 SINCE 1995	 ACCREDITED Management Systems Certification Body MSCB-156
This certificate remains the property of KTR-CC This certificate is accredited by IAS which is IAF MLA joined corporation. If not received further management audit within a given period, this certificate will be canceled.	

# Appendix 5, Certificate of Medical Device GMP

원본대조필

인정번호(No.) :KTR-AABA-7907

의료기기 제조 및 품질관리 기준 적합인정서  
(Certificate of GMP)

■ 업체명/허가번호(Company name of Applicant / License No.)  
(주)중헌제약  
Joonghun Pharmaceutical Co.,Ltd.

■ 대표자 (Representative)  
윤석준 ( Sukjoon Yoon )

■ 업체 소재지 (Company address of Applicant)  
충청북도 청주시 흥덕구 옥산면 옥산산단로 28 , 28101  
28 Oksansandanro, Chungjoosi, Choongchungbookdo 28101, Korea

■ 제조소명 (Name of Manufacturer)  
제조사 : (주)중헌제약(Joonghun Pharmaceutical Co.,Ltd.)

■ 제조소 소재지 (Address of Manufacturer)  
제조사 : 충청북도 청주시 흥덕구 옥산면 옥산산단로 28, 28101  
28 Oksansandanro, Chungjoosi, Choongchungbookdo 28101, Korea

■ 품목군 (Category)  
인체조직 또는 기능 대체품(Substitute of Human Tissue or Body Function)  
의료기기 제조 및 품질관리기준에 적합함을 인정합니다.  
(We hereby certify that the above manufacturer complies with Korea  
Good Manufacturing Practices of Medical Devices for the product group listed above)

발행일자(Date of Issue) 2020. 03. 27  
유효기간(Date of Expiration) : 2023. 03. 26

KTR

대전지방식품의약품안전청장  
DAEJEON REGIONAL FOOD AND DRUG ADMINISTRATION

한국화학융합시험연구원장  
KOREA TESTING & RESEARCH INSTITUTE


Headquarters  
Manufacturing Facility  
Tel.  
E-mail  
Website

15, Gukhoe-daero 62-gil, Yeongdeungpo-gu, Seoul, Republic of Korea, 07236  
28 Oksansandan-ro, Heungdeok-gu, Cheongju-si, Chungcheongbuk-do, Republic of Korea, 28101  
+82 2 2138 8114  
COVID-19TF@joonghun.com  
www.joonghun.com

Joonghun  
Pharmaceutical



## Appendix 6, Certificate of Pharmaceutical GMP



번호(No.) : MFDS-6-F-2631-1-2019-65

### 의약품 제조 및 품질관리기준 적합판정서 (Certificate of GMP Compliance of a Manufacturer)

- 제조소의 명칭/업허가번호(Name of Manufacturer/License No.)  
(주)충헌제약/2631
- 제조소의 소재지(Address of Manufacturer)  
충청북도 청주시 흥덕구 옥산면 옥산산단로 28 공장동
- 제형 또는 제조 방법  
(Dosage forms of Product(s) of Manufacturing operations : 붙임 참고)

위 제조소는 「약사법」 제38조, 「의약품 등의 안전에 관한 규칙」 제48조 및 의약품실사상호협력기구(PIC/S) 규정에 의한 의약품 제조 및 품질관리기준에 적합함을 판정합니다.

(We hereby certify that the above manufacturer complies with Good Manufacturing Practices of Pharmaceutical Products(s) according to the Korea Pharmaceutical Affairs Act and PIC/S GMP guides)

발급일자(Date of Issue) : 2019.11.06

유효기간(Date of Expiration) : 2022.04.16

비고 : 실사종료일 2019.04.17, 2019.06.19(주사제)

대전지방식품의약품안전청장

